

UPDATED: URGENT Field Safety Notice

Trilogy Evo, Trilogy Evo O2, Trilogy Evo Universal, Trilogy EV300
Non-pneumatic Nebulizer Use, Flow Sensor Nebulized Aerosol Deposition, and Obstruction Alarm

April 24, 2026

<To: Name / Title / Customer Name
Street Address
City, State, Zip Code
<modify title block format as needed>

This document contains important information for the continued safe and proper use of your equipment

Please review the following information with all members of your staff who need to be aware of the contents of this communication. It is important to understand the implications of this communication.

Please retain this letter for your records.

Dear Customer,

Philips Respironics is issuing this update to Field Safety Notice (FSN) 2026-CC-SRC-002, which was previously distributed regarding Trilogy Evo Platform ventilators (Trilogy Evo, Trilogy Evo O2, Trilogy Universal and Trilogy EV300 ventilators).

Since the initial communication, Philips Respironics has received feedback and questions from customers regarding the restriction on the use of non-pneumatic nebulizers (e.g., vibrating mesh nebulizers) as described in prior FSN 2026-CC-SRC-002 under reference bullet 2025-CC-SRC-020. Based on additional internal review and clinical evaluation, Philips Respironics has evaluated the risk versus benefit of continued use of the non-pneumatic nebulizers (e.g., vibrating mesh nebulizers). Based on this evaluation, Philips Respironics is revising its guidance to allow for the use of non-pneumatic nebulizers as long as the updated instructions for nebulizer placement and understanding of the warning contained in this Field Safety Notice are followed. Adherence to the revised guidance contained in this FSN will support continuity of care where alternative options may be limited or not appropriate in all care settings.


Therefore, this update provides clarification and revised guidance regarding the use of non-pneumatic nebulizers with Trilogy Evo Platform ventilators.

Specifically:

- Similar to pneumatic nebulizers, non-pneumatic nebulizers can be used as long as they are placed in accordance with the guidance and instructions on the proper location for a nebulizer on the various types of patient circuits. See “Figure 1: Acceptable In-Line Nebulizer Placement”

- **WARNING:** Ventilator accuracy can be adversely affected by changes to the inspired gas caused by the use of both pneumatic and non-pneumatic nebulizers.
- For customers who have updated 1.05.15.00 release, the User Interface only refers to pneumatic nebulizers – these instructions also apply to non-pneumatic nebulizers.

Figure 1: Acceptable In-Line Nebulizer Placement



This update **does not change the previously communicated information** regarding:

- Nebulized aerosol deposition on the internal flow sensor (FSN 2024-CC-SRC-013 – Flow Sensor Nebulized Aerosol Deposition), and
- Obstruction Alarm timing (IPN 2024-CC-SRC-002 – Obstruction Alarm).

Please review this letter carefully, as it contains **updated guidance regarding nebulizer use**. The specific nebulizer model should be selected and provided by the healthcare provider according to medication requirements. Philips Respironics cannot make clinical recommendations. Any questions regarding patient management should be directed to the appropriate clinical professionals. A response form is included and must be completed and returned as instructed to confirm receipt and understanding of this update.

1. What the problems are and under what circumstances they can occur

Use of non-pneumatic (e.g., Vibrating Mesh) Nebulizers (C&R 2025-CC-SRC-020)

Following the issuance of our Field Safety Notice regarding the use of Non-Pneumatic (e.g., Vibrating Mesh) Nebulizers with the Trilogy Evo device, we have become aware of the patient and operational impact associated with the restriction. At the time of the initial FSN issuance, the restriction on non-pneumatic (vibrating mesh) nebulizers was based on observed impacts to ventilator performance, including potential changes in delivered therapy accuracy due to aerosolized medication affecting the inspired gas.

Following evaluation of the risk versus benefit and customer feedback on real-world clinical, patient, and operational impact, Philips Respironics is revising its guidance to allow for the use of non-pneumatic nebulizers as long as the updated instructions and warning contained in this Field Safety Notice are followed.

2. Hazards and harms associated with nebulizer use (pneumatic or non-pneumatic)

Exercise caution when using a nebulizer (pneumatic or non-pneumatic) with the Trilogy Evo, Trilogy Evo O2, Trilogy Evo Universal and Trilogy EV300. Nebulizers may impact ventilator performance depending on the type of nebulizer, circuit configuration and therapy mode.

When a nebulizer (pneumatic or non-pneumatic) is used inline with a passive circuit during AC-VC, SIMV or AVAPS modes, clinicians should be aware of possible discrepancies in tidal volume delivery and monitoring. In adult active PAP and adult dual limb circuits, non-pneumatic nebulizers have not demonstrated impact on ventilator performance or displayed parameters based on available testing. In contrast, pneumatic nebulizers may affect ventilator performance across all circuit types.

During nebulization, volume delivery accuracy in AC-VC, SIMV, or AVAPS modes may be impacted by gas added by pneumatic nebulizers and, in certain configurations, changes in gas properties associated with non-pneumatic nebulizers. As a result, patient volume and flow monitoring displayed on the device may not reflect actual patient values and triggering sensitivity may be affected. In addition, aerosolized particles may impact circuit components including external flow sensor (EFS) (active flow and dual limb circuits) and passive exhalation valve (single limb circuits), mask or mask elbow. These discrepancies in delivery and monitoring can lead to potential harms such as hypoventilation, patient discomfort, and delayed clinical response. Any reduction in ventilator accuracy is expected to resolve once treatment is complete.

3. Affected products and how to identify them

All Trilogy Evo, Trilogy Evo O2, Trilogy Evo Universal and Trilogy EV300 models could be impacted by this letter.

4. Actions that should be taken by the customer or user in order to prevent risks for patients or users

- Follow the updated instructions contained within this notice regarding the use of non-pneumatic nebulizers
- Communicate this Urgent Field Safety Notice to anyone in your organization who interacts with a Trilogy Evo Platform ventilator, including clinicians, therapists, nurses, caregivers, and patients. Ensure they are aware of the changes regarding the use of nebulizers with these ventilators and that non-pneumatic nebulizers may be used.
- If a Trilogy Evo Platform ventilator previously in your possession has been transferred, ensure this information is communicated to them.
- Complete the URGENT Field Safety Notice Response Form included later in this letter. This is required to verify that you have received and understand the content of this notice.
- Distributors should identify customer list and where appropriate, distribute this Field Safety Notification and all relevant appendices to physicians, clinicians, patients and/or users.

5. Actions planned by Philips Respironics to correct the problem

Philips Respironics will be issuing an update to the User Manual to include this guidance on the use of Non-Pneumatic (e.g., Vibrating Mesh) Nebulizers. In addition, Philips Respironics will be releasing a software update to reflect this guidance on the user interface to clarify that it applies to both, non-pneumatic and pneumatic nebulizers. The updates will be included in a future release to minimize user impact through repeated device and literature updates.


If you need any further information or support concerning this issue, please contact your local Philips Respironics Customer Service at:

- For DMEs and Homecare customers: 1-800-345-6443, press option 2 for repair and then option 1
- For Hospital customers: 1 (800) 722-9377, press option 2

This notice has been reported to the appropriate Regulatory Agencies.

Philips regrets any inconvenience caused by this problem.

Sincerely,



Tracie Capozzio
Head of Quality Therapy Platforms,
Sleep and Respiratory Care

This addendum applies to the Trilogy Evo platform (Trilogy Evo, Trilogy Evo O2, and Trilogy EV300) and the Clinical and Caregiver manuals. The following sections are supplemental to or replace information within those manuals. The device operator is responsible for reading and understanding this addendum before use.



Device Usage

Contraindications

If the patient has any of the following conditions, consult the patient's health care professional before using noninvasive ventilation:

- An inability to maintain a patent airway or adequately clear secretions
- At risk to aspirate gastric contents
- Acute sinusitis or otitis media
- Epistaxis, causing pulmonary aspiration of blood
- Hypotension

The AVAPS-AE therapy mode is contraindicated for invasive use and patients less than 10 kg.

The AVAPS feature is contraindicated for patients less than 10 kg.

Users

- (Trilogy Evo and Trilogy Evo O2 only) Properly trained family, and non-professional caregivers.
- Clinical Users of the ventilator include: clinicians (respiratory and non-respiratory) and physicians.
- Service and Maintenance Users of the ventilator include: hospital and service provider technicians.

Warnings



Warning: Use only the cleaning methods outlined in this addendum. Philips is unable to verify the safety or performance of any device if ozone or other unapproved cleaning and disinfection methods are used.

- You must use a Philips Respironics-approved main flow bacteria filter on the patient gas outlet port to prevent patient or device contamination. Filters not approved by Philips Respironics may degrade system performance. For a list of accessories, see your device's accessories guide.
- Move the device away from any potential sources of electromagnetic interference (EMI), including MRI equipment, medical imaging systems, security systems, appliances, wireless communications equipment (such as cellular phones), computers, and televisions.



- Before placing a patient on the ventilator, perform a clinical assessment. Considerations should include:
 - Choosing alarm settings
 - Whether alternative ventilation equipment is required
 - Whether alternative monitors are required, such as Vte monitoring for Active PAP circuit, pulse oximeter, respiratory monitor with alarm, SpO₂, FiO₂, EtCO₂, and pulse rate.
- When using the ActivePAP circuit, CO₂ monitoring is required to measure exhaled carbon dioxide. Be certain that any CO₂ monitoring equipment used complies with ISO 80601-2-55.
- Never use an extension cord with the system. Always operate the device using a properly grounded AC power outlet. If you are unsure whether a power outlet is properly grounded, contact an electrician for assistance.

Notes

- Any serious incident that has occurred in relation to this device should be reported to Philips and the competent authority of the Member State in which the user and/or patient is established.


MRI Safety Information



The ventilation system is MR Unsafe. Keep it outside the MRI scan room (Zone IV). It represents a projectile hazard.

Symbols Glossary

See <http://www.symbols.philips.com> for a description of the symbols used on this device and its packaging.

Symbol	Definition
	Eject media

How to Contact Philips Respironics

Trilogy Evo, Trilogy Evo O2

Should you experience trouble with this equipment or require assistance, please contact your home care provider. If you need to contact Philips Respironics directly, contact customer service at +1-724-387-4000, or go to www.respironics.com to find your local customer service contact information.

Trilogy EV300

Should you experience trouble with this equipment or require assistance, contact Philips Respironics Customer Service. Visit www.healthcare.philips.com to find your local customer service contact information.



Device Overview

The following section is supplemental to the "About" chapter in your device manual.

Parts of the User Interface

Status Bar

Note: Features and icons may differ based on software version, product type, and/or market availability.



Use the status bar to monitor device status and the availability of manual therapeutic actions.

1	Deliver 100% Oxygen	5	Full access indicator
2	100% Oxygen timer	6	Bluetooth
3	Automatic algorithm restart	7	Bluetooth data transfer
4	Nebulizer Note: This icon will be displayed if the Nebulizer option is selected/set to Yes under the Device Options menu. Note: For more information, refer to the "Adding Nebulization" subsection, located within the "Device Setup (Startup)" section, of this addendum.	8	Alarm silence
		9	Power sources and their status
		10	Device Actions menu
		11	System time

Therapy Modes & Controls

Therapy Control Settings

Setting Name	Setting Range/Increment
Max Pressure	6-50 cmH ₂ O, 1 cmH ₂ O increments (as stated in device manual) Note: A maximum pressure setting of less than 40 cmH ₂ O is recommended for pediatric patients.
PC Min/Max	A/C-PC AVAPS: 0-40 cmH ₂ O Increments: 1 cmH ₂ O
PSV Min/Max	PSV AVAPS and AVAPS-AE: 0-40 cmH ₂ O Increments: 1 cmH ₂ O

Device Setup

The following section is supplemental to the "Device Setup" chapter in your device manual.

Connecting a Circuit

Leak Compensation

When using a Passive circuit, the ventilator provides leak compensation for intentional leak in the patient circuit and leaks that occur at the patient interface, such as cuff leak or mask leak. There is no leak compensation when the Active PAP, Active Flow, or Dual Limb circuits are used. Refer to the "Device Setup" chapter in your Clinical manual for more information about leak compensation.



Adding Nebulization



Warning: Do not use a non-pneumatic nebulizer, such as vibrating mesh, with this device. Doing so could degrade device performance, leading to harms such as hypoventilation, low oxygen saturation, and/or dyspnea.




Warning: When using a pneumatic nebulizer, ensure the nebulizer flow rate does not exceed 10 SLPM. Doing so could degrade device performance, leading to harms such as over-inflation of the lung, low oxygen saturation, hypoventilation, and/or delay/absence of therapy.

Correct Nebulizer Placement




Warning: Do not connect a pneumatic nebulizer to the ventilator breathing system (VBS) in any position other than as recommended in these instructions. Doing so could degrade device performance, leading to harms such as over-inflation of the lung, low oxygen saturation, hypoventilation, and/or delay/absence of therapy.

Passive





Active

(external flow sensor shown)



Dual Limb

(before or after wye connector)



If using pneumatic nebulizer treatments with the ventilator:

- The nebulizer should be installed **only** at the end of the patient connection port. For Dual Limb circuits, this could also be at the end of the patient's inspiratory limb, specifically at the input to the wye connector.
- Ensure the nebulizer is placed upright and below the patient tubing.
- Do not place a Heat Moisture Exchanger (HME) between the nebulizer and the patient interface.
- Prior to initiating nebulizer therapy, tap the **Nebulizer** button in the status bar to confirm correct nebulizer placement.

Note: The icon does not initiate nebulization therapy. Instead, it displays a prompt regarding correct nebulizer placement for the selected circuit type.



Nebulizer/Circuit Compatibility



Warning: Do **not** use a pneumatic nebulizer with Mouthpiece Ventilation (MPV) circuits. Doing so could degrade device performance, leading to harms such as over-inflation of the lung, low oxygen saturation, hypoventilation, and/or the delay/absence of therapy.

All **heated** circuits are compatible with pneumatic nebulizer use. However, some **non-heated** circuits and tidal volume combinations are not compatible, as described in the table below.



Warning: Do **not** use non-heated circuits in combination with pneumatic nebulizers that are listed as "NOT compatible" in the table below. Doing so could degrade device performance, leading to harms such as over-inflation of the lung, low oxygen saturation, hypoventilation, and/or delay/absence of therapy.

Pneumatic Nebulizer/Non-Heated Circuit Compatibility		
	Pediatric (14-18 mm)	Adult (20-22 mm)
<i>Passive</i>	Only compatible with tidal volume <300 mL	Only compatible with tidal volume <700 mL
<i>Active PAP</i>	NOT compatible	Compatible
<i>Active Flow</i>	NOT compatible	Compatible
<i>Dual Limb</i>	Compatible	Compatible


Adding Oxygen

Replacing Oxygen Fitting

The Trilogy Evo Oxygen Blending Module (OBM) is supplied from the factory with a DISS interface to allow connection to an oxygen source. Not all countries or facilities use this interface.

To replace the installed fitting:

1. Remove screws from OBM retaining plate.
2. Remove retaining plate.
3. Firmly pull oxygen fitting straight out of pocket.
4. Insert desired fitting into pocket with flat sides aligned. Fitting flange should be flush with face of pocket after installation.
5. Replace retaining plate with countersinks facing outward.
6. Replace screws to secure retaining plate.



Device Operation

New Patient

1. In the Home window, tap the **New Patient** button. This button clears all existing patient data. For more information about this button, refer to the following sections in your applicable clinical manual: "Setting Up a New Patient" (Trilogy Evo) or "Entering New Patient Information" (Trilogy EV300); and "Preparing the Device for Use by a Different Patient" (both).
2. Follow the prompts on the screen, and select and/or enter information as appropriate.
3. Acknowledge the reminder to ensure both a bacterial/viral filter and a particulate filter are installed.
4. Edit the prescription settings according to the procedure in the "About Prescriptions" (Trilogy Evo) or "About Prescriptions/Settings" (Trilogy EV300) section of your clinical manual.

Alarm & System Messages

Ventilator Inoperative

Priority	High
What to do	<p>Assess the patient and provide an alternate method of ventilation, then contact customer service.</p> <p>If the display screen is not operational, follow these steps:</p> <ol style="list-style-type: none"> 1. Hold the Alarm Silence and the On/Off (Standby) buttons for 5 seconds. 2. After 5 seconds, an audible tone will sound. Release both buttons. 3. Press the On/Off (Standby) button again. The device will stop delivering therapy and shut down.

Obstruction

Priority	High
Why it occurs	<p>The ventilator detects an obstruction in the patient's inhalation path, exhalation path, or external flow sensor. The ventilator detects that the leak device is missing.</p> <p>Note: The Obstruction alarm will not be triggered if there is an obstruction after the exhalation port (passive circuit) or after the active exhalation valve (Active PAP circuit). To detect obstructions after the exhalation port or active exhalation valve, use alternative alarms such as: High Inspiratory Pressure, Circuit Disconnect, Low Tidal Volume, Low Minute Ventilation, Low Respiratory Rate, Leakage, Rebreathing, Apnea, or Low Peak Inspiratory Pressure (volume modes) alarms.</p>
What to do	<p>Assess the patient, then:</p> <ul style="list-style-type: none"> • Check the circuit. Is it kinked or pinched? • Check the bacterial/viral filter and/or HME. Is it blocked? • Is the leak device blocked or missing? • Is the external flow sensor blocked?
Device performance	The device activates the obstruction alarm, evaluates the alarm condition until resolution of the obstruction, and continues to attempt to deliver therapy.



Algorithm summary	<p>The obstruction alarm can be triggered in one of the following ways:</p> <ul style="list-style-type: none"> • Two consecutive breaths have volumes less than a predetermined threshold based on circuit type and therapy settings. • Two consecutive breaths have average flows ≤ 1 LPM. • The average flow over 4 breath cycles is ≤ 1 LPM. • The average flow over 60 seconds is ≤ 1 LPM. <p>Additionally, when the external flow sensor is available:</p> <ul style="list-style-type: none"> • Two consecutive breaths have an absolute value of the average inspiratory flow ≤ 1 LPM and a flow span ≤ 3 LPM. • Flow remains ≤ 1 LPM for 60 seconds.
-------------------	---

Circuit Disconnect

Priority	High
What to do	<p>Assess the patient, then:</p> <ul style="list-style-type: none"> • Is the circuit connected to the patient? • Is the circuit connected to the ventilator? • Does a large unplanned leak exist? • Is the tracheostomy or endotracheal tube occluded?

High Tidal Volume

Priority	Medium
Why it occurs	<p><i>Passive, active flow, or dual limb circuit types:</i> The estimated exhaled tidal volume is more than or equal to the alarm setting for a number of breaths depending on the therapy mode:</p> <ul style="list-style-type: none"> • Three consecutive breaths: A/C-PC, CPAP, PSV, S/T, SIMV-VC, SIMV-PC • Six consecutive breaths: A/C-VC <p>When AVAPS or AVAPS-AE mode is enabled, this alarm occurs when the Tidal Volume is more than or equal to the alarm threshold for one minute and Set Inspiratory Pressure is less than or equal to (Set Inspiratory Pressure min + 1 cmH₂O).</p> <p><i>Active PAP circuit type:</i> The delivered tidal volume is more than or equal to the alarm setting for three consecutive breaths.</p>



Low Tidal Volume

Priority	Medium
Why it occurs	<p><i>Passive, active flow, or dual limb circuit types:</i> The estimated exhaled tidal volume is less than or equal to the low tidal volume alarm setting for a number of breaths depending on the therapy mode:</p> <ul style="list-style-type: none"> • Three consecutive breaths: A/C-PC, CPAP, PSV, S/T, SIMV-VC, SIMV-PC • Six consecutive breaths: A/C-VC <p>When AVAPS or AVAPS-AE mode is enabled, this alarm occurs when the Tidal Volume is less than or equal to the alarm setting for one minute and Set Inspiratory Pressure is greater than or equal to (Set Inspiratory Pressure max - 1 cmH₂O).</p> <p><i>Active PAP circuit type:</i> The delivered tidal volume is less than or equal to the low tidal volume setting.</p>

Inlet Filter(s) Blocked

Priority	Low
Why it occurs	The inlet filter(s) becomes blocked and delivered therapy is reduced.
What to do	<ul style="list-style-type: none"> • Is the air inlet(s) blocked? • Remove, rinse, and reinstall the air-inlet foam filter. See "Rinsing the Air-Inlet Foam Filter" in your device manual. • Discard and replace the particulate filter.
Device performance	The device continues to function.

Patient-related Alarm Availability by Therapy Mode

Alarm	A/C-PC	A/C-VC	MPV-PC*	MPV-VC*	CPAP	PSV	S/T	SIMV-PC	SIMV-VC	AVAPS-AE
Apnea (requires Backup Ventilation)					•	•		•	•	
High Inspiratory Pressure		•		•					•	
Low Inspiratory Pressure		•		•					•	
No Trigger			•	•						
High SpO ₂	•	•	•	•	•	•	•	•	•	•
Low SpO ₂	•	•	•	•	•	•	•	•	•	•
High Pulse Rate	•	•	•	•	•	•	•	•	•	•
Low Pulse Rate	•	•	•	•	•	•	•	•	•	•

*MPV-PC and MPV-VC available for Trilogy Evo and Trilogy Evo O2 only.



Testing Circuit Obstruction Alarms

The following section is supplemental to the “Testing Circuit Obstruction Alarms” section in your device manual. For ventilator-dependent patients, do not rely on any single alarm to detect when a circuit is obstructed. One or more of the following alarms may indicate an obstructed circuit.

- Obstruction
- High Inspiratory Pressure
- Circuit Disconnected
- Low Tidal Volume
- Low Minute Ventilation
- Low Respiratory Rate
- Low Peak Inspiratory Pressure alarms (settable for volume modes)
- Leakage alarm (Active PAP circuit only)
- Rebreathing Detected

To test that these alarms detect a circuit obstruction, follow these steps.

1. **For a Passive circuit:** Remove the leak device, and block the end of the circuit.

For an Active PAP circuit: Remove the active exhalation valve (AEV), and block the end of the circuit.

For an Active Flow or a Dual Limb circuit: Block the end of the circuit.


2. Start ventilation.
3. Confirm that one or more of the alarms listed above activate.
4. Stop ventilation.

Disposal

Dispose of this device in accordance with local collections and recycling regulations. For more information, visit www.philips.com/recycling.

Accessories and Replacement Parts

Prior to use, inspect all cords, filters, and circuits for damage or wear. Discard and replace if necessary.



Bacterial/Viral Filter



Particulate Filter




Air-Inlet Foam Filter



Bacterial/Viral Filter

The use of a bacterial/viral filter is always required on the patient gas outlet port.



Particulate Filter

The use of the particulate filter is required to protect the device from dirt and dust. For instructions on changing the particulate filter, see Replacing the Particulate Filter in the Service and Maintenance section of this addendum.

- In the hospital setting, replace at least monthly.
- When used outside of the hospital setting, replace at least monthly and between patients.

Note: The filter is printed with the **manufacturing** date (FAB [MM/YYYY]); the filter must be used within 5 years of that date.

External Flow Sensor (Active Flow & Dual Limb Circuits)

Note: For optimal performance, ensure sensor is always in an upright position (i.e., cap and cable are on top).

Humidification Accessories

If you are using a humidifier with the ventilator, refer to the instructions provided with the humidifier.

Detachable Battery

- **Warning:** Do not incinerate, disassemble, short circuit or expose to temperatures above 60°C (140°F). Please refer to the instructions for use.
- **Warning:** Do not disassemble, strike, crush or throw into fire. In the event of a serious bout, do not continue to use. Do not place in high temperature environment. Do not use the battery after immersion in water!

FiO₂ Sensor

Calibration

Prior to calibration, the device should be turned on and/or plugged into AC power for a minimum of 15 minutes. Be sure the air-inlet foam filter and the particulate filter are installed in the device. Refer to your device manual for additional information.


Installation

Required: T20 Torx screwdriver




To install a new FiO₂ sensor:

1. Using a T20 Torx screwdriver, unscrew the FiO₂ sensor access panel and set aside. Note that the panel is also the FiO₂ sensor driver, used in step 3.




2. Place the sensor in the FiO₂ socket.



3. Place the FiO₂ driver on the sensor and twist clockwise to secure the sensor in the receptacle.



4. Push the connector onto the pins on the sensor. Be sure you have aligned the pins with the receptacles on the connector.



5. Replace the sensor access panel and secure with the screws.

Removal

Required: T20 Torx screwdriver and a slender, pointed tool such as a scratch awl

To remove the FiO₂ sensor:

1. Using a T20 Torx screwdriver, unscrew the FiO₂ sensor access panel and set aside.
2. Use a slender, pointed tool such as a scratch awl to pull the tab lock away from the connector to the unlock position.
3. While holding the tab in the unlock position, pull the white connector straight out. Do not use the wires to pull the connector.
4. Place the FiO₂ driver (sensor access panel) on the sensor and twist counterclockwise to remove the sensor from the device.

USB Flash Drive

The USB flash drive (removable storage) included with your device can be used to import and export data.



Cleaning & Disinfection

Note: Any damage caused by unapproved ozone cleaning and disinfection methods or other unapproved cleaning and disinfection methods will not be covered by the Philips Limited Warranty.

Note: If a home-based device is returned and assigned to a new patient, clean and disinfect the device exterior according to the procedures in this section.

The following section replaces the information presented in the Cleaning and Disinfection chapter in your device manual.

Warnings

- Because the device is intended for multi-patient use, be sure you follow the cleaning and disinfection instructions in this addendum.
- You must use a Philips Respironics-approved main flow bacteria filter on the patient gas outlet port to prevent patient or device contamination. Filters not approved by Philips Respironics may degrade system performance. For a list of accessories, see your device's accessories guide.
- This product is suitable for use with other patients if a bacteria filter was in use at all times and exchanged between patients.
- There are no approved air path disinfection methods. Use only the cleaning methods outlined in this addendum. Philips is unable to verify the safety or performance of any device if ozone or other unapproved cleaning and disinfection methods are used.

Exterior Cleaning and Disinfection

Warning: To avoid electric shock, do not remove the enclosure cover. Only service personnel should remove the enclosure. After cleaning and disinfecting, be sure the device is completely dry before reattaching accessories and connectors and before reconnecting it to a power source. To avoid electrical shock, always unplug the power cord from the wall outlet before cleaning the ventilator. If the device has been exposed to rain or dampness, dry the device (including the area around the power cord connection) with the power cord disconnected from the device before applying AC power.

Caution: Do not immerse the device or allow liquids into any of the controls or inside the enclosure as the device may be damaged. If this occurs, contact your equipment provider for assistance. Use only the cleaning agents and methods described in this section to clean and disinfect the device.

Cleaning the Exterior

Frequency: Clean the device's exterior surface weekly and between patients.

Requirements:

- Lint-free cloth
- Soft-bristle brush
- Liquid dishwashing detergent solution: 1 teaspoon (5 mL) of liquid dishwashing detergent (such as Dawn Ultra Dishwashing Liquid®) per gallon (3.7 L) of potable water

To clean the exterior, follow these steps:

1. Turn the device off and disconnect it from the power source.
2. Detach all accessories and connectors.
3. Use a lint-free cloth dampened (not dripping) with a liquid dishwashing detergent solution to clean the exterior of the enclosure.
4. Use a soft-bristle brush in the areas around the screen, buttons, and any other areas where soil may be difficult to remove. Ensure you remove all visible soil.
5. Use a lint-free cloth dampened (not dripping) with potable water to remove all detergent residue.



6. Use a lint-free cloth to dry the enclosure.
7. Inspect the device for cleanliness.
8. Repeat the cleaning steps until the surfaces are visibly clean.
9. Inspect the device for damage after cleaning. If any parts are damaged, contact customer service.

Disinfecting the Exterior

Frequency: Disinfect the exterior surface weekly and between patients.

Prerequisite: Before disinfecting the exterior, be sure you have cleaned the device as instructed in the previous section, "Cleaning the Exterior."

Isopropyl Alcohol

Requirement: 70% isopropyl alcohol, lint-free cloth

To disinfect with alcohol, follow these steps:

1. Use a lint-free cloth dampened with alcohol to wipe the alcohol onto the exterior, thoroughly wetting the surfaces.
2. Keep wet 10 minutes.
3. Allow to air dry.
4. Inspect the device for damage after disinfection. If any parts are damaged, contact customer service.

Chlorine Bleach

Requirement: Household chlorine bleach containing 8.25% sodium hypochlorite, lint-free cloth

To disinfect with bleach, follow these steps:

1. Combine 10 parts potable water to 1-part bleach.
2. Use a lint-free cloth dampened with the bleach solution to wipe the bleach solution onto the exterior, thoroughly wetting the surfaces.
3. Keep wet 10 minutes.
4. Allow to air dry.
5. Inspect the device for damage after disinfection. If any parts are damaged, contact customer service.

Air Path Cross Contamination

Warning: You must use a Philips Respironics-approved main flow bacteria filter on the patient gas outlet port to prevent patient or device contamination. Filters not approved by Philips Respironics may degrade system performance.

To manage multiple patients, you can select either the use of a bacteria filter or replacement of the gas pathway components in between patient use.

This product is suitable for multi-patient use when used with a bacteria filter to reduce cross contamination between patients. Ensure the bacteria filter is replaced between patients per the instructions provided with the bacteria filter. If a previous patient has not used a bacteria filter, the next patient must use a bacteria filter, or the gas pathway components must be replaced. For more information about exchanging the gas pathway components, contact Philips Respironics Customer Service at 1-724-387-4000 or visit www.respironics.com for your local customer service number.

Note: *The bacteria filter should be disposed of in accordance with the filter instructions for use and between patients.*




Service and Maintenance

Replacing the Air-Inlet Foam Filter

The air-inlet foam filter is the gray foam located on the back panel. It is required to protect the device from dirt and dust. Replace if it becomes damaged.

In the clinical environment, replace monthly and between patients. In the home environment, replace every six months and between patients. Only use Philips Respironics-supplied filters. Dispose according to local regulations. Ventilation can continue while you are replacing the filter.

Note: The air-inlet foam filter may be replaced by the user, caregiver, clinician, or service personnel.



To replace the disposable inlet filter:

1. Be sure you have a replacement filter nearby.
2. Pinch the filter and pull it out of the filter cover.
3. Insert the clean replacement filter into the filter cover. Be sure it is positioned securely.


Replacing the Particulate Filter

The use of the particulate filter is required to protect the device from dirt and dust. Ventilation can continue while you are replacing the filter.


- In the hospital setting, replace at least monthly.
- When used outside of the hospital setting, replace at least monthly and between patients.

Note: The particulate filter may be replaced by the user, caregiver, clinician, or service personnel.


Twist the filter cover counterclockwise a quarter of a turn, and then pull straight out to remove.



Twist the filter counterclockwise a quarter of a turn, and then pull straight out to remove.



Place a new filter onto the bayonet mount then twist the filter clockwise a quarter of a turn while pressing in to secure.



Replace the filter cover and turn clockwise to secure.

Technical Data

Measured and Displayed Patient Parameters

Pulse rate with pulse oximeter accessory	Nonin Accessories: 18 to 321 beats per minute in increments of 1 beat per minute Masimo Accessories: 25 to 240 beats per minute in increments of 1 beat per minute See "Pulse Oximeter."
EPAP	Expiratory positive airway pressure. Displayed parameter in AVAPS-AE mode.

Environmental

Operating	<ul style="list-style-type: none"> • Temperature: 0°C to 40°C, 32°F to 104°F • Relative humidity: 5% RH to 90% RH, non-condensing • Atmospheric pressure: 62 kPa to 106 kPa • Altitude: Approximately -384 m to 3954 m (-1,253 ft to 12,971 ft) • Battery charging temperature: 5°C to 40°C, 32°F to 104°F
-----------	---



Wireless

NFC	Note: This function is currently not available.
Wi-Fi	Note: This function is currently not available.

Pulse Oximeter

Displayed pulse rate range:	Nonin Accessories: 18 to 321 beats per minute in increments of 1 beat per minute Masimo Accessories: 25 to 240 beats per minute in increments of 1 beat per minute
Data averaging	Nonin Accessories: 4 beat average, updated every second Masimo Accessories: 8 second average, updated every second

Regulatory Information

Wireless Communication

Bluetooth Core Specification version 4.2

Wireless

Hereby, Respironics Inc. declares that this class 1 radio equipment is in compliance with Directive 2014/53/EU. The full text of the EU declaration of conformity is available at the following internet address:
<http://incenter.medical.philipscom/PMSPublic>.

This product meets RF exposure requirements when it is positioned with a separation distance of at least 20 cm away from the body.

Note: The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Philips Respironics is under license.



Limited Warranty

The following information replaces the information presented in the Warranty chapter in the Trilogy Evo, Trilogy Evo O2, and Trilogy EV300 manuals.

Respironics, Inc., a Philips company ("Philips Respironics") provides this non-transferable, limited warranty for Trilogy Evo, Trilogy Evo O2, and Trilogy EV300 ("Product") to the customer who originally purchased the Product directly from Philips Respironics.

What this warranty covers: Philips Respironics warrants each new Product will be free from defects in materials and workmanship and will perform in accordance with the Product specifications under normal and proper use and maintenance in accordance with applicable instructions, subject to the exclusions below.

How long does this warranty last: Two years from the date of shipment by Philips Respironics to the original purchaser (Trilogy Evo, Trilogy Evo O2), or one year from the longer of the date of shipment to the purchaser or date of setup by purchaser for the end user (Trilogy EV300), except:

- a. The warranty period for the internal battery included with the Product is 90 days from the date of shipment to the original purchaser.
- b. All other accessories and replacement parts are not covered under this warranty.

What this warranty does not cover: This warranty does not apply to any software included with the Product as the software warranty is included in the software license. This warranty does not cover damage or injury whether to the Products, personal property, or persons caused by accident, misuse, abuse, Acts of God, water ingress, unapproved ozone cleaning and disinfection methods, other unapproved cleaning and disinfection methods, repair or alteration by anyone other than Philips Respironics or its authorized service center, failure to operate in accordance with the terms of the operating manual and instructions, lack of reasonable care, the discontinuance of a network (e.g., 2G, 3G, etc.) by a carrier (e.g., ATT, Verizon, etc.), or other defects not related to material or workmanship. This warranty is not transferable. If Philips Respironics finds that a Product returned for service or the issue raised is not covered under this limited warranty, Philips Respironics may charge an evaluation fee and return shipping.

What Philips Respironics will do: If a Product fails to conform to the warranties set forth above during the applicable warranty period, Philips Respironics will repair or replace the Product or refund the original purchase price, in Philips Respironics sole discretion. Philips Respironics may use new or remanufactured assemblies, components, and parts in repair and new or recertified refurbished devices for replacement. The balance of the original warranty period will apply to any Product or component of a Product repaired or replaced under this warranty.

Warranty Disclaimer; Limitation of Liability: EXCEPT AS SET FORTH IN THIS LIMITED WARRANTY, PHILIPS RESPIRONICS MAKES NO WARRANTIES, EXPRESSED OR IMPLIED, STATUTORY OR OTHERWISE, REGARDING THE PRODUCT OR ITS QUALITY OR PERFORMANCE. PHILIPS RESPIRONICS SPECIFICALLY DISCLAIMS THE IMPLIED WARRANTY OF MERCHANTABILITY AND THE IMPLIED WARRANTY OF FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL PHILIPS RESPIRONICS MAXIMUM LIABILITY UNDER THESE WARRANTIES EXCEED THE ORIGINAL PURCHASE PRICE OR WILL PHILIPS RESPIRONICS BE LIABLE FOR ANY ECONOMIC LOSS, LOSS OF PROFITS, OVERHEAD, OR SPECIAL, INCIDENTAL, OR CONSEQUENTIAL DAMAGES. Repair, replacement, or return of purchase price by Philips Respironics is the original purchaser's sole and exclusive remedy under this warranty.

This warranty gives you specific legal rights, and you may also have other rights that vary from state to state. Some states do not allow the exclusion or limitation of incidental or consequential damages, so the above exclusion and limitations may not apply to you.

How to get warranty support: Patients contact your local authorized Philips Respironics dealer and dealers contact Respironics, Inc. at:

1001 Murry Ridge Lane
Murrysville, Pennsylvania 15668-8550
+1-724-387-4000



Note: For Australian and New Zealand customers, this warranty replaces the warranty included with the user manual.

1. The following statement is provided to a customer who is a consumer under the Australian Consumer Law: Our goods come with guarantees that cannot be excluded under the Australian Consumer Law. You are entitled to a replacement or refund for a major failure and for compensation for any other reasonably foreseeable loss or damage. You are also entitled to have the goods repaired or replaced if the goods fail to be of acceptable quality and the failure does not amount to a major failure. This guarantee applies in addition to other rights and remedies available to the consumer under law.
2. The following statement is provided to a customer who is a consumer under the Consumer Guarantees Act 1993 (NZ): For consumers under the Consumer Guarantees Act 1993 (NZ) who purchase the goods for personal, domestic or household use: Our goods come with guarantees that cannot be excluded under the Consumer Guarantees Act 1993 (NZ).
3. Philips warrants that the products shall be free from defects of workmanship and materials, and will perform in accordance with the product specifications for a period of two (2) years from the date of shipment by Philips Respironics to the original purchaser (Trilogy Evo, Trilogy Evo O2), or one (1) year from the longer of the date of shipment to the purchaser or date of setup by purchaser for the end user (Trilogy EV300). This Warranty covers the replacement or repair at the option of Philips, of any product that has a manufacturing or material defect that is not the result of normal wear and tear, or a natural characteristic of the material used. This Warranty is not transferable and does not cover products used for commercial purposes, and it does not apply to any consumable items (including but not limited to filters, masks, tubes and humidifier chambers).
4. The customer is responsible for returning the product to an authorised Philips Homecare Provider, and for collecting the product from the Homecare Provider after repair or replacement, at its own cost. Philips is responsible only for the freight cost of transporting the product between the Homecare Provider and Philips. Philips reserves the right to charge an evaluation and postage fee for any returned product where no problem is found following evaluation.
5. This Warranty does not cover:
 - products purchased outside of Australia and New Zealand;
 - any damage caused as a result of misuse or abuse, modification, tampering with or alteration of the product, pest infestation, liquid egress into the product, or unapproved ozone cleaning and disinfection methods;
 - contamination due to cigarette, pipe, cigar or other smoke;
 - failure to follow manufacturer’s instruction for use as per the user manual;
 - defects that are a consequence of repairs to a product made or attempted by a service provider other than one approved by Philips;
 - products that have been subjected to incorrect electrical supply or inconsistent electrical supply or used with inappropriate accessories.
6. This Warranty is not transferrable in the event of any resale or transfer of products.
7. To make a claim under this Warranty, contact your Homecare Provider or Philips:

AUSTRALIA
Philips Electronics Australia Limited
65 Epping Road, North Ryde
NSW 2113 Australia
Tel: 1300 766 488
Email: repairs-src@philips-easyconnect.com

NEW ZEALAND
Philips New Zealand Commercial Ltd
Level 3, 123 Carlton Gore Road
New Market Auckland 1023 New Zealand
Tel: 0800 251 400
Email: repairs-src@philips-easyconnect.com



Appendix A -- uSpO2 Masimo SET Oximetry Cable, USB

Intended Use

The Masimo SET uSpO2 Pulse Oximetry Cable is indicated for the continuous noninvasive monitoring of functional oxygen saturation of arterial hemoglobin (SpO2) and pulse rate (measured by an SpO2 sensor). The Masimo SET uSpO2 Pulse Oximetry Cable is indicated for use with adult, pediatric and neonatal patients during both motion and no motion conditions, and for patients who are well or poorly perfused in hospitals and hospital-type facilities.

Compatibility

These accessories are compatible with the Trilogy Evo and Trilogy EV300 series of ventilators when used in a hospital environment. They are also compatible with LNCS Adhesive and LNCS Reusable sensors. Refer to the sensor's instructions for accuracy.

Warnings

- The pulse oximeter is to be operated by, or under the supervision of, qualified personnel only. The manual, accessories, directions for use, all precautionary information, and specifications should be read before use.
- Do not start or operate the pulse oximeter unless the setup was verified to be correct.
- Explosion hazard: Do not use the pulse oximeter in the presence of flammable anesthetics or other flammable substance in combination with air, oxygen-enriched environments, or nitrous oxide.
- To ensure safety, avoid stacking multiple devices or placing anything on the device during operation.
- To protect from electric shock, always remove the sensor and completely disconnect the pulse oximeter before bathing the patient.
- If any measurement seems questionable, first check the patient's vital signs by alternate means and then check the pulse oximeter for proper functioning
- The pulse oximeter is not an apnea monitor.
- The pulse oximeter may be used during electrocautery, but this may affect the accuracy or availability of the parameters and measurements.
- The pulse oximeter should not be used for arrhythmia analysis.

Cautions

- Electrical shock and flammability hazard: Before cleaning, always turn off the device and disconnect from any power source.
- Do not place the pulse oximeter on electrical equipment that may affect the device, preventing it from working properly.
- To ensure that alarm limits are appropriate for the patient being monitored, check the limits each time the pulse oximeter is used.
- Variation in measurements may be profound and may be affected by sampling technique as well as the patient's physiological conditions. Any results exhibiting inconsistency with the patient's clinical status should be repeated and/or supplemented with additional test data. Blood samples should be analyzed by laboratory instruments prior to clinical decision making to completely understand the patient's condition.
- Do not submerge the pulse oximeter in any cleaning solution or attempt to sterilize by autoclave, irradiation, steam, gas, ethylene oxide or any other method. This will seriously damage the pulse oximeter.
- To minimize radio interference, other electrical equipment that emits radio frequency transmissions should not be in close proximity to the pulse oximeter.

Notes

- A functional tester cannot be used to assess the accuracy of the pulse oximeter.



- Do not loop the patient cabling into a tight coil or wrap around the device, as this can damage the patient cabling.
- Additional information specific to the Masimo sensors compatible with the pulse oximeter, including information about parameter/measurement performance during motion and low perfusion, may be found in the sensor's directions for use (DFU).

Instructions for Use

For additional information, warnings, and instructions to use the sensor or cable, see that accessory's instructions for use. To attach the sensor or cable to the ventilator, see the ventilator's instructions for use.

Environmental Specifications

For environmental specifications, see the accessory's instructions for use.

Cleaning

For cleaning instructions, see the accessory's instructions for use.

Disposal

Dispose according to local regulations. If you need help, contact Philips Respironics.

Specifications

Displayed Oxygen Saturation Range (SpO ₂):	0 to 100% with a resolution of 1%
Displayed pulse rate range:	25 to 240 beats per minute with a resolution of 1
SpO ₂ and pulse rate accuracy:	See the sensor instructions.
Data update period:	Every second
Data averaging:	8 second average, updated every second
Radiant Power at 50 mA pulsed:	≤15 mW

The Masimo SET technology with Masimo sensors has been validated for no motion accuracy in human blood studies on healthy adult male and female volunteers with light to dark skin pigmentation in induced hypoxia studies in the range of 70-100% SpO₂ against a laboratory CO-Oximeter and ECG monitor. This variation equals ±1 standard deviation. Plus or minus one standard deviation encompasses 68% of the population.

The Masimo SET technology with Masimo sensors has been validated for motion accuracy in human blood studies on healthy adult male and female volunteers with light to dark skin pigmentation in induced hypoxia studies while performing rubbing and tapping motions, at 2 to 4 Hz at an amplitude of 1 to 2 cm and a non-repetitive motion between 1 to 5 Hz at an amplitude of 2 to 3 cm in induced hypoxia studies in the range of 70-100% SpO₂ against a laboratory CO-Oximeter and ECG monitor. This variation equals ±1 standard deviation, which encompasses 68% of the population.

The Masimo SET technology has been validated for low perfusion accuracy in bench top testing against a Biotek Index 2™ simulator and Masimo's simulator with signal strengths of greater than 0.02% and transmission of greater than 5% for saturations ranging from 70 to 100%. This variation equals ±1 standard deviation. Plus or minus one standard deviation encompasses 68% of the population.

The Masimo SET Technology with Masimo Neo sensors has been validated for neonatal motion accuracy in human blood studies on healthy adult male and female volunteers with light to dark skin pigmentation in induced hypoxia studies while performing rubbing and tapping motions, at 2 to 4 Hz at an amplitude of 1 to 2 cm and a non-repetitive motion between 1 to 5 Hz at an amplitude of 2 to 3 cm in induced hypoxia studies in the range of 70-100% SpO₂ against a laboratory CO-Oximeter and ECG monitor. This variation equals ±1 standard deviation. Plus or minus one standard deviation encompasses 68% of the population. 1% has been added to the results to account for the effects of fetal hemoglobin present in neonates.



The Masimo SET technology with Masimo sensors has been validated for pulse rate accuracy for the range of 25 -240 bpm in bench top testing against a Biotek Index 2™ simulator. This variation equals ± 1 standard deviation. Plus or minus one standard deviation encompasses 68% of the population.

See sensor directions for use (DFU) for complete application information. Unless otherwise indicated, reposition reusable sensors at least every 4 hours and adhesive sensors at least every 8 hours.

Sensor accuracy specified when used with Masimo technology using a Masimo patient cable for LNOP sensors, RD SET sensors, the LNCS sensors, or the M-LNCS sensors. Numbers represent Arms (RMS error compared to the reference). Because pulse oximeter measurements are statistically distributed, only about two-thirds of the measurements can be expected to fall within a range of \pm Arms compared to the reference value. Unless otherwise noted, SpO2 accuracy is specified from 70% to 100%. Pulse Rate accuracy is specified from 25 to 240 BPM.

Masimo oximetry sensors use light emitting diodes (LEDs) which operate at the following wavelengths (except TC-I, TF-I, and TFA-1 sensor types as specified in their respective columns below.)

LED	General	TC-I	TF-I	TFA-1
	Nominal Wavelength	Nominal Wavelength	Nominal Wavelength	Nominal Wavelength
Red	660 nm	653 nm	660 nm	653 nm
Infrared	905 nm	880 nm	880 nm	905 nm

EMC Information

Note: The EMISSIONS characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.

End-User License Agreement

THIS DOCUMENT IS A LEGAL AGREEMENT BETWEEN YOU ("PURCHASER") AND PHILIPS. IF YOU DO NOT AGREE TO THE TERMS OF THIS AGREEMENT, PROMPTLY RETURN THE ENTIRE PACKAGE, INCLUDING ALL ACCESSORIES, IN THEIR ORIGINAL PACKAGE, WITH YOUR SALES RECEIPT TO PHILIPS FOR A FULL REFUND.

1. Grant of License. In consideration of payment of the license fee, which is part of the price paid for this product, Philips grants to Purchaser a nonexclusive, nontransferable license, without right to sublicense, to use the copy of the incorporated software/firmware, and documentation in connection with Purchaser's use of the Masimo Products for their labeled purpose. Philips reserves all rights not expressly granted to Purchaser.
2. Ownership of Software/Firmware. Title to, ownership of, and all rights and interests in, any Masimo software and/or firmware and the documentation, and all copies thereof, remain at all times vested in Masimo Corporation, licensor to Philips, and they do not pass to Purchaser.
3. Assignment. Purchaser shall not assign or transfer this License, in whole or in part, by operation of law or otherwise, without Philips' prior written consent; any attempt without such consent, to assign any rights, duties or obligations arising hereunder shall be void.
4. Copy Restrictions. The software/firmware, mask works, circuit board layouts, and accompanying written materials are copyrighted. Unauthorized copying of the software, including software that has been modified, merged, or included with other software, or other written materials is expressly forbidden. You may be held legally responsible for any copyright infringement that is cause or incurred by your failure to abide by the terms of this license. Nothing in this license provides any rights beyond those provided by 17 U.S.C. §117.
5. Use Restriction. As the Purchaser, you may physically transfer the products from one location to another provided that the software/firmware is not copied. You may not electronically transfer the software/firmware from the products to any other device. You may not disclose, publish, translate, release distribute copies of, modify, adapt, translate,



reverse engineer, decompile, disassemble, or create derivative works based on the Masimo Product, the software/firmware, or the written materials without the prior written consent of Masimo. Masimo Sensors that are designated for single use are licensed under Masimo patents for use on a single patient only, and are not sold. Masimo reusable sensors are licensed under Masimo patents for use on more than a single patient, but are not sold. Possession of this device does not convey any express or implied license to use the device with unauthorized sensors or cables that would, alone, or in combination with this device, fall within the scope of one or more of the patents relating to this device. There is no license, implied or otherwise, that would allow use of single use Masimo Sensors beyond their intended single use. After use of single use Masimo Sensors, there is no further license granted by Masimo to use the sensors and they must be discarded.

6. Non-Authorized Accessories. Masimo technology is designed to operate together with Masimo cables, sensors and accessories as an integrated system. When any component of the system is compromised, erroneous measurements can occur. Accordingly, the use of unauthorized sensors or accessories, such as third-party reprocessed or copycat sensors can yield unreliable results when used with a Masimo device. The performance of Masimo technology is not validated when used with any unauthorized sensor or accessory.

7. Transfer Restrictions. The software/firmware is licensed to the Purchaser, and may not be transferred to anyone, except other end-users, without the prior written consent of Philips. In no event may you transfer, assign, rent, lease, sell, or otherwise dispose of the software/firmware or the products on a temporary basis.

8. Beneficiary. Masimo Corporation is a Beneficiary of this Agreement and has the right to enforce its provisions.

9. U.S. Government Rights. If you are acquiring software (including the related documentation) on behalf of any part of the United State Government, the following provisions apply: the software is deemed to be "commercial software" and "commercial computer software documentation," respectively pursuant to DFAR Section 227.7202 FAR 12.212, as applicable. Any use, modification, reproduction, release, performance, display or disclosure of the software (including the related documentation) by the U.S. Government or any of its agencies shall be governed solely by the terms of this Agreement and shall be prohibited except to the extent expressly permitted by the terms of this agreement.

10. Patents. The Masimo technology in this product is covered under one or more patents as set forth at www.masimo.com/patents.htm

